

Genesis Touch® 3.4

Product Overview

With its robust feature set, the Genesis Touch® 3.4 leverages an Android™ tablet for an intuitive, consumer-based interface, connected to a full range of LifeStream™ applications, services and devices.

What's New in Genesis Touch® 3.4:

- Compatibility with Welch Allyn *Bluetooth*® peripherals, including a weight scale and two non-invasive blood pressure monitors.
 - Welch Allyn RPM-Scale10
 - Welch Allyn H-BP100SBP (1700 Series)
 - Welch Allyn RPM-BP100 (1500 Series)
- Faster *Bluetooth*® peripheral pairing process.

The Genesis Touch® 3.4 offers:

- Wireless *Bluetooth*® connectivity to a growing variety of peripheral medical devices.
- Integrated video to improve the quality of care and enable dynamic patient interaction.
- Educational content to encourage patient compliance and healthy behaviors.
- A user-friendly platform that easily integrates with LifeStream™.
- Canadian-French and Spanish language support for patients and care providers.

Integrated Video Improves Quality of Care and Enables Patient Interaction

- Enable patients to easily join a video visit with their care providers.
- Multiple care providers and family members may interact with patients during a video visit.
- Increase educational effectiveness with optional call recording and playback.
- Improve patient interaction while providing lower-cost access to patients in hard-to-reach areas.
- Optional 4G video-over-cellular data plans promote patient mobility and stable video clarity.
- Leverage a flexible, cost-effective licensing structure to alleviate the need for patient licensing.



This personal health device helps care providers foster meaningful patient relationships, by enabling simple video-based communication in addition to traditional telemonitoring.

A patient education feature encourages healthier choices and compliance with a personalized care plan.

Support for peripheral devices meets the needs of patients with different risk factors and specialized monitoring requirements.

Genesis Touch® 3.4 Technical Specifications

An Emphasis on Education

- Replay programmed questions on demand, not only during scheduled monitoring time.
- Questions are self-guided and address topics like diet, diabetes, cardiac care, medications & more.
- Encourage patients to become more involved with managing their own health.
- Combine with video visits for face-to-face education opportunities.

Mobility Enhances Patient Comfort, Care and Convenience

- Access patients in hard-to-reach rural areas or patients with limited Internet access.
- Synchronize with a patient's lifestyle and routine with optimal in-home monitor placement.
- Reach a mobile population with Wi-Fi and 4G/3G cellular connectivity.

A Secure Data Stream: From the Home to the Hospital

- Class two personal health device promotes stability, safety, and confidentiality.
- Dedicated, secure health device, allows care providers to enable HIPAA compliant operations.
- Stable, tested and controlled interface protects against viruses and patient tampering.
- Integrated connectivity status menu to troubleshoot & download remote configurations instantly.

A Common Platform Enables Growth and Accessibility

- Leverage the consumer-based Android™ platform.
- Interoperable with a growing variety of interfaces and applications.
- Support for manual entry of vital signs like PT/INR, peak flow and blood glucose encourages patients to continue using a familiar device.

SPECIFICATIONS

Leverage a Samsung Galaxy Tab® E platform, powered by Android™ 6.0.1, 4G cellular LTE™ and Wi-Fi.

A VARIETY OF PERIPHERAL MEDICAL DEVICES

Reduce installation times with wireless *Bluetooth*® peripherals and no physical wiring in residence. Compatibility with the following medical devices helps reach more patients with a variety of risk factors.

Scales:

Honeywell HomMed® Scale
A&D UC355PBT-Ci Scale
Accuro HRS305 Handrail Scale*
Welch Allyn RPM-Scale10

Blood Pressure Monitor

A&D UA767PBT-Ci Monitor
Welch Allyn H-BP100SBP (1700 Series)*
Welch Allyn RPM-BP100 (1500 Series)*

Pulse Oximeters:

Nonin Onyx® II 9560
CONTEC CMS50E
CONTEC CMS50EW*

Glucometers:

Bayer® CONTOUR® NEXT EZ
Bayer® CONTOUR®
Johnson & Johnson OneTouch Ultra®
Johnson & Johnson OneTouch Ultra®2
Johnson & Johnson OneTouch UltraMini®
Entra Health MyGlucoHealth® MGH-BT1

For more information

www.honeywelllifecare.com

Honeywell Life Care Solutions

3400 Intertech Drive
Suite 200
Brookfield, WI 53045
888-353-5404
www.honeywell.com

HomMed®, LifeStream™ and Genesis Touch® are either trademarks or registered trademarks of Honeywell HomMed®, LLC.

Android™ is a trademark of Google Inc.; The Samsung Galaxy Tab® platform is a trademark of SAMSUNG; Welch Allyn, Connex, SureTemp, FlexiPort, and SureBP are registered trademarks of Welch Allyn. Nonin Onyx® is a registered trademark of Nonin Medical, Inc.; Wi-Fi is a registered certification mark of the Wi-Fi Alliance; Bayer® and the Bayer Cross® are trademarks of Bayer. CONTOUR® is a trademark of Ascensia Diabetes Care US, Inc.; OneTouch and Ultra® are registered trademarks of Johnson & Johnson; MyGlucoHealth® is a registered trademark of Entra Health Systems; LTE™ is a registered trademark of ETSI. Fuze is a trademark of the FuzeBox Software Corporation.

All other product names, trademarks, logos and brands cited herein, are property of their respective owners. All company, product and service names used in this document are for identification purposes only.

*NCI Technology, Inc., Accuro Scale; Welch Allyn H-BP100SBP; Welch Allyn RPM-BP100; CONTEC CMS50EW devices are not available in Canada